




DETMOLD D-95 Surgical Respirator Mask fitting overview



Securing respirators with the D95 + Band + FitClip combination

The D95+ Surgical Respirator is designed to achieve a comfortable and adequate seal by joining the mask's two bands with the FitClip™. Available in two sizes, the FitClip™ allows flexibility in the adequate fitting of respirators.

The FitClip™ offers a secure hold point on one end, and an open end for ease of attaching and removing the second mask band.

The Detmold Medical range of D95+ Surgical Respirators have the Standard Fitclip™ already attached, and where a longer band FitClip is required, the XL FitClip™ can be ordered separately.

An inadequate seal between the Respirator and the face can result in contaminated air entering the mask, please refer to your appropriate Respirator Fit Test Procedure to determine the required combination of mask and FitClip™.

INSTRUCTIONS FOR USE

Prior to donning, ensure hand hygiene is undertaken

Applying the mask



Removing the mask safely



- 1.** Move the FitClip to the band on your preferred side for donning.
- 2.** Open the mask, flattening the nose piece.

- 3.** Lean chin into the mask and lift bands toward the crown of the head, securing onto the FitClip.
- 4.** Adjust the bands for a secure and comfortable fit

- 5.** Using both hands, firmly mould the nosepiece over bridge of the nose.
- 6.** Pinch the nosepiece (angled over the nostril) and drop the jaw gently allowing the mask to fit to the natural shape of your face.

- 7.** Look left to right, then up and down, ensuring there is no movement in the mask at the nosepiece.
- 8.** Inhale and exhale deeply checking for escape of air. If air is escaping, refit the mask.

- 1.** Without touching the mask, follow bands to the FitClip and gently unhook, keeping tension on the bands.
- 2.** Lean forward maintaining control of the band tension, whilst slowly removing the mask away from the face. Ensure the mask does not flick forward or back onto the face.
- 3.** Discard safely

There should be no gaps between the face and mask, providing an adequate seal. Once the mask is applied it should not be touched or reapplied. Appropriate hand hygiene should always be used